



Meatballs CL, Ready to Cook

INGREDIENTS:

100 kg	Pork	S II
25 kg	Drinking water for swelling	
12 kg	Meatball Fix CL	Art. No. 536201
0.5 kg	Nubaron CL/AF	Art. No. 436500

PROCESSING:

Allow the Meatball Fix CL to swell in the drinking water for 10 minutes. Combine the pork with Nubaron CL/AF and mince to a 5 mm grain size. Then mix the minced meat with the pre-swelled Meatball Fix CL until well combined and allow to set. Then portion and shape.

PREPARATION:

Fry on both sides in a frying pan or on the grill.

ATTENTION:

In its raw state, this product falls under Regulation (EC) 853/2004!

NON-BINDING DECLARATION RECOMMENDATION:

Pork, drinking water, potato fiber, rice flour, spices (contains MUSTARD), table salt, WHOLE EGG POWDER, starch, seasoning, dextrose, yeast extract, spice extracts

Recipe ID: 180

* All compound foods and ingredients used must be labelled with those of the respective manufacturer in addition to our information.